

Mouse Casp14/Caspase-14 ELISA Kit

Catalog No: #EK19949

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Casp14/Caspase-14 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Caspase-14
Target Name	Caspase-14
Other Names	CASP-14,Mini-ICE,MICE
Accession No.	O89094
Uniprot	O89094
GeneID	12365;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39.2pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Casp14

Note: This product is for in vitro research use only