

Mouse Amotl1/Angiotensin-like protein 1 ELISA Kit

Catalog No: #EK19950

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Amotl1/Angiotensin-like protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Angiotensin-like protein 1
Target Name	Angiotensin-like protein 1
Other Names	junction-enriched and -associated protein, JEAP
Accession No.	Q9D4H4
Uniprot	Q9D4H4
GeneID	75723;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.3pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Amotl1

Note: This product is for in vitro research use only