

Human AMOTL2/Angiotensin-like protein 2 ELISA Kit

Catalog No: #EK19951

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human AMOTL2/Angiotensin-like protein 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Angiotensin-like protein 2
Target Name	Angiotensin-like protein 2
Other Names	Leman coiled-coil protein, LCCP, KIAA0989
Accession No.	Q9Y2J4
Uniprot	Q9Y2J4
GeneID	51421;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: AMOTL2

Note: This product is for in vitro research use only