

Mouse Amotl2/Angiotensin-like protein 2 ELISA Kit

Catalog No: #EK19952

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Amotl2/Angiotensin-like protein 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Angiotensin-like protein 2
Target Name	Angiotensin-like protein 2
Accession No.	Q8K371
Uniprot	Q8K371
GeneID	56332;

Application Details

Assay Type: Sandwich Detect Range: 0.625-40ng/mL Sensitivity: 0.327ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Amotl2

Note: This product is for in vitro research use only