

Mouse Ccr2/C-C chemokine receptor type 2 ELISA Kit

Catalog No: #EK19953

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ccr2/C-C chemokine receptor type 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse C-C chemokine receptor type 2
Target Name	C-C chemokine receptor type 2
Other Names	C-C CKR-2, JE/FIC receptor, MCP-1 receptor, CD192, Cmkbr2
Accession No.	P51683
Uniprot	P51683
GeneID	12772;

Application Details

Assay Type: Sandwich Detect Range: 0.625-40ng/mL Sensitivity: 0.315ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Ccr2

Note: This product is for in vitro research use only