

Rat Ccr2/C-C chemokine receptor type 2 ELISA Kit

Catalog No: #EK19954

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Ccr2/C-C chemokine receptor type 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat C-C chemokine receptor type 2
Target Name	C-C chemokine receptor type 2
Other Names	C-C CKR-2,CD192,Cmkbr2
Accession No.	O55193
Uniprot	O55193
GeneID	60463;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ccr2

Note: This product is for in vitro research use only