

Human C1GALT1/Glycoprotein-N-acetylgalactosamine 3-beta-galactosyltransferase 1 ELISA Kit

Catalog No: #EK19955

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human C1GALT1/Glycoprotein-N-acetylgalactosamine 3-beta-galactosyltransferase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Glycoprotein-N-acetylgalactosamine 3-beta-galactosyltransferase 1
Target Name	Glycoprotein-N-acetylgalactosamine 3-beta-galactosyltransferase 1
Other Names	B3Gal-T8,Core 1 O-glycan T-synthase,Core 1 UDP-galactose:N-acetylgalactosamine-alpha-R beta 1,3-galactosyltransferase 1,Beta-1,3-galactosyltransferase,Core 1 beta1,3-galactosyltransferase 1,C1GalT1
Accession No.	Q9NS00
Uniprot	Q9NS00
GeneID	56913;

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.052ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:C1GALT1

Note: This product is for in vitro research use only