

Mouse C1galt1/Glycoprotein-N-acetylgalactosamine 3-beta-galactosyltransferase 1 ELISA Kit



Catalog No: #EK19956

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse C1galt1/Glycoprotein-N-acetylgalactosamine 3-beta-galactosyltransferase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Glycoprotein-N-acetylgalactosamine 3-beta-galactosyltransferase 1
Target Name	Glycoprotein-N-acetylgalactosamine 3-beta-galactosyltransferase 1
Other Names	Core 1 O-glycan T-synthase,T-syn,Core 1 UDP-galactose:N-acetylgalactosamine-alpha-R beta 1,3-galactosyltransferase 1,Core 1 beta1,3-galactosyltransferase 1,C1GalT1,Pit1
Accession No.	Q9JJ06
Uniprot	Q9JJ06
GeneID	94192;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.091ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:C1galt1

Note: This product is for in vitro research use only