

Human AHI1/Jouberin ELISA Kit

Catalog No: #EK19958

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human AHI1/Jouberin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Jouberin
Target Name	Jouberin
Other Names	Abelson helper integration site 1 protein homolog,AHI-1
Accession No.	Q8N157
Uniprot	Q8N157
GeneID	54806;

Application Details

Assay Type: Sandwich Detect Range:31.2-2000pg/mL Sensitivity:15.62pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:AHI1

Note: This product is for in vitro research use only