

Rat Atg12/Ubiquitin-like protein ATG12 ELISA Kit

Catalog No: #EK19959

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Atg12/Ubiquitin-like protein ATG12 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Ubiquitin-like protein ATG12
Target Name	Ubiquitin-like protein ATG12
Other Names	Autophagy-related protein 12,APG12-like,Apg12l
Accession No.	Q2TBJ5
Uniprot	Q2TBJ5
GeneID	361321;

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/ml Sensitivity:0.407ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Atg12

Note: This product is for in vitro research use only