

Human AXL/Tyrosine-protein kinase receptor UFO ELISA Kit



Catalog No: #EK19960

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human AXL/Tyrosine-protein kinase receptor UFO ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Tyrosine-protein kinase receptor UFO
Target Name	Tyrosine-protein kinase receptor UFO
Other Names	AXL oncogene,UFO
Accession No.	P30530
Uniprot	P30530
GeneID	558;

Application Details

Assay Type: Sandwich Detect Range:31.2-2000pg/mL Sensitivity:13.7pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:AXL

Note: This product is for in vitro research use only