

Mouse Axl/Tyrosine-protein kinase receptor UFO ELISA Kit

Catalog No: #EK19961

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Axl/Tyrosine-protein kinase receptor UFO ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Tyrosine-protein kinase receptor UFO
Target Name	Tyrosine-protein kinase receptor UFO
Other Names	Adhesion-related kinase,Ark,Ufo
Accession No.	Q00993
Uniprot	Q00993
GeneID	26362;

Application Details

Assay Type: SandwichDetect Range:0.625-40ng/mL Sensitivity:0.329ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Axl

Note: This product is for in vitro research use only