

# Mouse B4galt5/Beta-1,4-galactosyltransferase 5 ELISA Kit



Catalog No: #EK19962

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse B4galt5/Beta-1,4-galactosyltransferase 5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Beta-1,4-galactosyltransferase 5
Target Name	Beta-1,4-galactosyltransferase 5
Other Names	Beta-1,4-GalTase 5, Beta-1,4-GalT II, UDP-Gal:beta-GlcNAc beta-1,4-galactosyltransferase 5, UDP-galactose:beta-N-acetylglucosamine beta-1,4-galactosyltransferase 5, Bgt-5
Accession No.	Q9JMK0
Uniprot	Q9JMK0
GeneID	56336;

## Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/ML Sensitivity: 7.8pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: B4galt5

Note: This product is for in vitro research use only