

Mouse Cblb/E3 ubiquitin-protein ligase CBL-B ELISA Kit



Catalog No: #EK19963

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cblb/E3 ubiquitin-protein ligase CBL-B ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse E3 ubiquitin-protein ligase CBL-B
Target Name	E3 ubiquitin-protein ligase CBL-B
Other Names	Casitas B-lineage lymphoma proto-oncogene b,RING-type E3 ubiquitin transferase CBL-B,SH3-binding protein CBL-B,Signal transduction protein CBL-B
Accession No.	Q3TTA7
Uniprot	Q3TTA7
GeneID	208650;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:33 pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cblb

Note: This product is for in vitro research use only