

## Rat Cd2ap/CD2-associated protein ELISA Kit

Catalog No: #EK19965

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Cd2ap/CD2-associated protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat CD2-associated protein
Target Name	CD2-associated protein
Accession No.	F1LRS8
Uniprot	F1LRS8

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.172ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cd2ap

Note: This product is for in vitro research use only