

## Human CD5L/CD5 antigen-like ELISA Kit

Catalog No: #EK19966

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CD5L/CD5 antigen-like ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human CD5 antigen-like
Target Name	CD5 antigen-like
Other Names	Apoptosis inhibitor expressed by macrophages,hAIM,CT-2,IgM-associated peptide,SP-alpha,API6,UNQ203/PRO229
Accession No.	O43866
Uniprot	O43866
GeneID	922;

## Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CD5L

Note: This product is for in vitro research use only