

# Human CTCF/Transcriptional repressor CTCF ELISA Kit



Catalog No: #EK19968

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CTCF/Transcriptional repressor CTCF ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Transcriptional repressor CTCF
Target Name	Transcriptional repressor CTCF
Other Names	11-zinc finger protein,CCCTC-binding factor,CTCF paralog
Accession No.	P49711
Uniprot	P49711
GeneID	10664;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.164ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CTCF

Note: This product is for in vitro research use only