

Rat Ctcf/Transcriptional repressor CTCF ELISA Kit

Catalog No: #EK19969

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Ctcf/Transcriptional repressor CTCF ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Transcriptional repressor CTCF
Target Name	Transcriptional repressor CTCF
Other Names	11-zinc finger protein,CCCTC-binding factor,CTCF paralog
Accession No.	Q9R1D1
Uniprot	Q9R1D1
GeneID	83726;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.168ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ctcf

Note: This product is for in vitro research use only