

Human CTSC/Dipeptidyl peptidase 1 ELISA Kit

Catalog No: #EK19970

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Human CTSC/Dipeptidyl peptidase 1 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Hu |
| Specificity | Natural and recombinant Human Dipeptidyl peptidase 1 |
| Target Name | Dipeptidyl peptidase 1 |
| Other Names | Cathepsin C, Cathepsin J, Dipeptidyl peptidase I, DPP-I, Dipeptidyl transferase, CPPI |
| Accession No. | P53634 |
| Uniprot | P53634 |
| GeneID | 1075; |

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CTSC

Note: This product is for in vitro research use only