

Mouse Ctsc/Dipeptidyl peptidase 1 ELISA Kit

Catalog No: #EK19971

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ctsc/Dipeptidyl peptidase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Dipeptidyl peptidase 1
Target Name	Dipeptidyl peptidase 1
Other Names	Cathepsin C, Cathepsin J, Dipeptidyl peptidase I, DPP-I, Dipeptidyl transferase
Accession No.	P97821
Uniprot	P97821
GeneID	13032;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.08ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Ctsc

Note: This product is for in vitro research use only