

Rat Ctsc/Dipeptidyl peptidase 1 ELISA Kit

Catalog No: #EK19972

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Ctsc/Dipeptidyl peptidase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Dipeptidyl peptidase 1
Target Name	Dipeptidyl peptidase 1
Other Names	Cathepsin C,Cathepsin J,Dipeptidyl peptidase I,DPP-I,Dipeptidyl transferase
Accession No.	P80067
Uniprot	P80067
GeneID	25423;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/ml Sensitivity:0.069ng/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ctsc

Note: This product is for in vitro research use only