

Mouse Crcp/DNA-directed RNA polymerase III subunit
RPC9 ELISA Kit

Catalog No: #EK19973



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse Crcp/DNA-directed RNA polymerase III subunit RPC9 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse DNA-directed RNA polymerase III subunit RPC9
Target Name	DNA-directed RNA polymerase III subunit RPC9
Other Names	RNA polymerase III subunit C9,Calcitonin gene-related peptide-receptor component protein,CGRP-RCP
Accession No.	O35427
Uniprot	O35427
GeneID	12909;

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.067ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Crcp

Note: This product is for in vitro research use only