

Rat Crcp/DNA-directed RNA polymerase III subunit RPC9 ELISA Kit



Catalog No: #EK19974

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Crcp/DNA-directed RNA polymerase III subunit RPC9 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat DNA-directed RNA polymerase III subunit RPC9
Target Name	DNA-directed RNA polymerase III subunit RPC9
Other Names	RNA polymerase III subunit C9, Calcitonin gene-related peptide-receptor component protein, CGRP-RCP
Accession No.	Q8VHM6
Uniprot	Q8VHM6
GeneID	114205;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.134ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Crcp

Note: This product is for in vitro research use only