

# Human CEMP1/Cementoblastoma-derived protein 1 ELISA Kit



Catalog No: #EK19975

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CEMP1/Cementoblastoma-derived protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cementoblastoma-derived protein 1
Target Name	Cementoblastoma-derived protein 1
Other Names	Cementum protein 1,Cementum protein 23,CP-23,CP23
Accession No.	Q6PRD7
Uniprot	Q6PRD7
GeneID	752014;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.157ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CEMP1

Note: This product is for in vitro research use only