

# Mouse Crispld2/Cysteine-rich secretory protein LCCL domain-containing 2 ELISA Kit



Catalog No: #EK19978

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Crispld2/Cysteine-rich secretory protein LCCL domain-containing 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cysteine-rich secretory protein LCCL domain-containing 2
Target Name	Cysteine-rich secretory protein LCCL domain-containing 2
Other Names	Coffeecrisp
Accession No.	Q8BZQ2
Uniprot	Q8BZQ2
GeneID	78892;

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.159ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Crispld2

Note: This product is for in vitro research use only