

Human ANGPTL6/Angiopoietin-related protein 6 ELISA Kit



Catalog No: #EK19979

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ANGPTL6/Angiopoietin-related protein 6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Angiopoietin-related protein 6
Target Name	Angiopoietin-related protein 6
Other Names	Angiopoietin-like protein 6,Angiopoietin-related growth factor,Angiopoietin-related protein 5,AGF,ARP5,UNQ152/PRO178
Accession No.	Q8NI99
Uniprot	Q8NI99
GeneID	83854;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.077ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ANGPTL6

Note: This product is for in vitro research use only