

# Mouse Angptl6/Angiopoietin-related protein 6 ELISA Kit



Catalog No: #EK19980

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Angptl6/Angiopoietin-related protein 6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Angiopoietin-related protein 6
Target Name	Angiopoietin-related protein 6
Other Names	Angiopoietin-like protein 6,Angiopoietin-related growth factor,Agf
Accession No.	Q8R0Z6
Uniprot	Q8R0Z6
GeneID	70726;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.2pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Angptl6

Note: This product is for in vitro research use only