

Human AHNAK2/Protein AHNAK2 ELISA Kit

Catalog No: #EK19981



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Human AHNAK2/Protein AHNAK2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Protein AHNAK2
Target Name	Protein AHNAK2
Other Names	C14orf78,KIAA2019
Accession No.	Q8IVF2
Uniprot	Q8IVF2
GeneID	113146;

Application Details

Assay Type: SandwichDetect Range:78-5000pg/mL Sensitivity:22pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:AHNAK2

Note: This product is for in vitro research use only