

# Rat Cited4/Cbp/p300-interacting transactivator 4 ELISA Kit

Catalog No: #EK19982

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                                   |
|--------------------|-------------------------------------------------------------------|
| Product Name       | Rat Cited4/Cbp/p300-interacting transactivator 4 ELISA Kit        |
| Brief Description  | ELISA Kits                                                        |
| Applications       | ELISA                                                             |
| Species Reactivity | Rt                                                                |
| Specificity        | Natural and recombinant Rat Cbp/p300-interacting transactivator 4 |
| Target Name        | Cbp/p300-interacting transactivator 4                             |
| Other Names        | MSG1-related protein 2,MRG-2,Mrg2                                 |
| Accession No.      | Q99MA0                                                            |
| Uniprot            | Q99MA0                                                            |
| GeneID             | 114491;                                                           |

## Application Details

Assay Type: SandwichDetect Range:78-5000pg/mL Sensitivity:39pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:Cited4

Note: This product is for in vitro research use only