

# Mouse Adamts18/A disintegrin and metalloproteinase with thrombospondin motifs 18 ELISA Kit



Catalog No: #EK19983

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Adamts18/A disintegrin and metalloproteinase with thrombospondin motifs 18 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse A disintegrin and metalloproteinase with thrombospondin motifs 18
Target Name	A disintegrin and metalloproteinase with thrombospondin motifs 18
Other Names	ADAM-TS 18
Accession No.	Q4VC17
Uniprot	Q4VC17
GeneID	208936;

## Application Details

Assay Type: Sandwich Detect Range: 0.625-40ng/mL Sensitivity: 0.303ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Adamts18

Note: This product is for in vitro research use only