

## Rat Cldn19/Claudin-19 ELISA Kit

Catalog No: #EK19985

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Cldn19/Claudin-19 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Claudin-19
Target Name	Claudin-19
Accession No.	Q5QT56
Uniprot	Q5QT56
GeneID	298487;

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.183ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cldn19

Note: This product is for in vitro research use only