

Human B4GALT6/Beta-1,4-galactosyltransferase 6 ELISA Kit

Catalog No: #EK19986

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human B4GALT6/Beta-1,4-galactosyltransferase 6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Beta-1,4-galactosyltransferase 6
Target Name	Beta-1,4-galactosyltransferase 6
Other Names	Beta-1,4-GalTase 6,UDP-Gal:beta-GlcNAc beta-1,4-galactosyltransferase 6,UDP-galactose:beta-N-acetylglucosamine beta-1,4-galactosyltransferase 6
Accession No.	Q9UBX8
Uniprot	Q9UBX8
GeneID	9331;

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.086ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:B4GALT6

Note: This product is for in vitro research use only