

Mouse Aox1/Aldehyde oxidase 1 ELISA Kit

Catalog No: #EK19987

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Aox1/Aldehyde oxidase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Aldehyde oxidase 1
Target Name	Aldehyde oxidase 1
Other Names	Azaheterocycle hydroxylase 1, Retinal oxidase, Ao, Ro
Accession No.	O54754
Uniprot	O54754
GeneID	11761;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10 ng/mL Sensitivity: 0.086 ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100 µl Assay Time: 3-5 h Detection wavelength: 450 nm gene_full_name: Aox1

Note: This product is for in vitro research use only