

Mouse AOX1_/Aldehyde oxidase ELISA Kit

Catalog No: #EK19988

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse AOX1_/Aldehyde oxidase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Aldehyde oxidase
Target Name	Aldehyde oxidase
Other Names	Retinal oxidase,Ao,Ro
Accession No.	O54754
Uniprot	O54754
GeneID	11761;

Application Details

Assay Type: Sandwich Detect Range:31.2-2000pg/mL Sensitivity:14.4pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:AOX1_

Note: This product is for in vitro research use only