

Rat Cnnm4/Metal transporter CNNM4 ELISA Kit

Catalog No: #EK19990



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rat Cnnm4/Metal transporter CNNM4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Metal transporter CNNM4
Target Name	Metal transporter CNNM4
Other Names	Ancient conserved domain-containing protein 4,Cyclin-M4,Acdp4
Accession No.	P0C588
Uniprot	P0C588
GeneID	363216;

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.053ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Cnnm4

Note: This product is for in vitro research use only