

Rat Cnm2/Metal transporter CNNM2 ELISA Kit

Catalog No: #EK19991

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cnm2/Metal transporter CNNM2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Metal transporter CNNM2
Target Name	Metal transporter CNNM2
Other Names	Ancient conserved domain-containing protein 2,Cyclin-M2,Acdp2
Accession No.	Q5U2P1
Uniprot	Q5U2P1
GeneID	294014;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.158 ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cnm2

Note: This product is for in vitro research use only