

## Human ATP1B4/Protein ATP1B4 ELISA Kit

Catalog No: #EK19992

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ATP1B4/Protein ATP1B4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Protein ATP1B4
Target Name	Protein ATP1B4
Other Names	X,K-ATPase subunit beta-m,X/potassium-transporting ATPase subunit beta-m
Accession No.	Q9UN42
Uniprot	Q9UN42
GeneID	23439;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:ATP1B4

Note: This product is for in vitro research use only