

Bovine APMAP/Adipocyte plasma membrane-associated protein ELISA Kit



Catalog No: #EK19994

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine APMAP/Adipocyte plasma membrane-associated protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Adipocyte plasma membrane-associated protein
Target Name	Adipocyte plasma membrane-associated protein
Accession No.	Q3T0E5
Uniprot	Q3T0E5
GeneID	535740;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:APMAP

Note: This product is for in vitro research use only