

Mouse Apmmap/Adipocyte plasma membrane-associated protein ELISA Kit

Catalog No: #EK19995

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Apmmap/Adipocyte plasma membrane-associated protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Adipocyte plasma membrane-associated protein
Target Name	Adipocyte plasma membrane-associated protein
Other Names	Protein DD16
Accession No.	Q9D7N9
Uniprot	Q9D7N9
GeneID	71881;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 34 pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Apmmap

Note: This product is for in vitro research use only