

Rat Apmmap/Adipocyte plasma membrane-associated protein ELISA Kit

Catalog No: #EK19996

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Apmmap/Adipocyte plasma membrane-associated protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Adipocyte plasma membrane-associated protein
Target Name	Adipocyte plasma membrane-associated protein
Other Names	Ab2-305
Accession No.	Q7TP48
Uniprot	Q7TP48
GeneID	366227;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/ml Sensitivity: 41pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Apmmap

Note: This product is for in vitro research use only