

Mouse Cd200/OX-2 membrane glycoprotein ELISA Kit

Catalog No: #EK19997

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cd200/OX-2 membrane glycoprotein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse OX-2 membrane glycoprotein
Target Name	OX-2 membrane glycoprotein
Other Names	MRC OX-2 antigen, CD200, Mox2
Accession No.	O54901
Uniprot	O54901

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/ml Sensitivity: 0.086ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cd200

Note: This product is for in vitro research use only