

Bovine CDK12/Cyclin-dependent kinase 12 ELISA Kit

Catalog No: #EK19998

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine CDK12/Cyclin-dependent kinase 12 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Cyclin-dependent kinase 12
Target Name	Cyclin-dependent kinase 12
Other Names	Cell division protein kinase 12
Accession No.	E1BB50
Uniprot	E1BB50

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.165ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CDK12

Note: This product is for in vitro research use only