

Rat Cdk12/Cyclin-dependent kinase 12 ELISA Kit

Catalog No: #EK20000

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cdk12/Cyclin-dependent kinase 12 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cyclin-dependent kinase 12
Target Name	Cyclin-dependent kinase 12
Other Names	Cdc2-related kinase, arginine/serine-rich, CrkRS, Cell division cycle 2-related protein kinase 7, CDC2-related protein kinase 7, Cell division protein kinase 12, Protein kinase for splicing component, Crk7, Crkrs, Pksc
Accession No.	Q3MJK5
Uniprot	Q3MJK5
GeneID	192350;

Application Details

Assay Type: Sandwich Detect Range: 78-5000 pg/mL Sensitivity: 42 pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100 µl Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cdk12

Note: This product is for in vitro research use only