

Mouse Cdh6/Cadherin-6 ELISA Kit

Catalog No: #EK20001

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cdh6/Cadherin-6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cadherin-6
Target Name	Cadherin-6
Other Names	Kidney cadherin, K-cadherin
Accession No.	P97326
Uniprot	P97326
GeneID	12563;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.087ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cdh6

Note: This product is for in vitro research use only