

Human CRTAC1/Cartilage acidic protein 1 ELISA Kit

Catalog No: #EK20002

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CRTAC1/Cartilage acidic protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cartilage acidic protein 1
Target Name	Cartilage acidic protein 1
Other Names	68 kDa chondrocyte-expressed protein, CEP-68, ASPIC, ASPIC1, CEP68
Accession No.	Q9NQ79
Uniprot	Q9NQ79
GeneID	55118;

Application Details

Assay Type: Sandwich Detect Range: 62.5-4000pg/mL Sensitivity: 31.25pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CRTAC1

Note: This product is for in vitro research use only