

## CCNK Antibody

Catalog No: #44708

Package Size: #44708-1 50ul #44708-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	CCNK Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	CCNK antibody detects endogenous levels of total CCNK
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CCNK
Target Name	CCNK
Other Names	CCNK antibody CCNK_HUMAN antibody CPR 4 antibody CPR4 antibody Cyclin K antibody Cyclin-K antibody CyclinK antibody MGC9113 antibody
Accession No.	Swiss-Prot#: O75909 NCBI Gene ID:8812
Uniprot	O75909
GeneID	8812;
Calculated MW	64 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: □ 1:500-1:2000

## Images

Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.



---

Note: This product is for in vitro research use only