

Human AEBP1/Adipocyte enhancer-binding protein 1 ELISA Kit



Catalog No: #EK20003

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human AEBP1/Adipocyte enhancer-binding protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Adipocyte enhancer-binding protein 1
Target Name	Adipocyte enhancer-binding protein 1
Other Names	AE-binding protein 1,Aortic carboxypeptidase-like protein,ACLP
Accession No.	Q8IUX7
Uniprot	Q8IUX7
GeneID	165;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:42.3pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:AEBP1

Note: This product is for in vitro research use only