



Mouse Arc/Activity-regulated cytoskeleton-associated protein ELISA Kit

Catalog No: #EK20004

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse Arc/Activity-regulated cytoskeleton-associated protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Activity-regulated cytoskeleton-associated protein
Target Name	Activity-regulated cytoskeleton-associated protein
Other Names	ARC/ARG3.1,mArc,Activity-regulated gene 3.1 protein homolog,Arg3.1
Accession No.	Q9WV31
Uniprot	Q9WV31
GeneID	11838;

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.171ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Arc

Note: This product is for in vitro research use only