

Mouse Adal/Adenosine deaminase-like protein ELISA Kit

Catalog No: #EK20006

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse Adal/Adenosine deaminase-like protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Adenosine deaminase-like protein
Target Name	Adenosine deaminase-like protein
Accession No.	Q80SY6
Uniprot	Q80SY6
GeneID	75894;

Application Details

Assay Type: SandwichDetect Range:3.12-200U/L Sensitivity:1.54U/L Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Adal

Note: This product is for in vitro research use only