

Mouse Actb/Actin, cytoplasmic 1 ELISA Kit

Catalog No: #EK20007

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Actb/Actin, cytoplasmic 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Actin, cytoplasmic 1
Target Name	Actin, cytoplasmic 1
Other Names	Beta-actin
Accession No.	P60710
Uniprot	P60710
GeneID	11461;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.12ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Actb

Note: This product is for in vitro research use only